


Advanced Light Microscopy Course, GSK Cancer Engineering

Alexandros Pertsinidis, PhD

Structural Biology Program


pertsina@mskcc.org

Imaging across spatial scales



What a microscope needs to do

- Magnify things
- Resolve points which are close together
- Collect as much light as possible (esp. for fluorescence)
- Do all of the above while introducing as little distortion as possible



Also, for imaging live biological samples

- **Speed is important**
- **Must not disturb biological process**

Image position and magnification depend on lens curvature (focal length) and on the physical distance from the object to the lens




Figure 4-3
Geometrical optics of a simple lens. The focal length f , focal point F , object-lens distance a , and lens-image distance b are indicated.




Figure 4-6
Object-image relationships. (a) $a < f$; (b) $a = f$; (c) $2f > a > f$; (d) $a = 2f$; (e) $a > 2f$.

All simple lenses have associated aberrations

- Still may encounter: *chromatic aberration* on cheap microscopes (prism effect--but can be reduced by using monochromatic light), *spherical aberration* when imaging deep into samples (e.g. embryos, even when the objective is "corrected"), *field curvature* when using bright lenses for fluorescence (but this is not a problem if you're imaging cells only in the center of the field)




Figure 4-8

Aberrations of a simple lens. (a) Chromatic aberration: Parallel incident rays of different wavelength are focused at different locations. (b) Spherical aberration: Incident rays parallel to the optic axis and reaching the center and the periphery of the lens are focused at different locations. (c) Coma: Off-axis rays passing through the center and periphery of the lens are focused at different locations. (d) Astigmatism: An off-axis aberration causes waves passing through the vertical and horizontal diameters to focus an object point as a streak. (e) Distortion and field curvature: The image plane is curved and not planar. So-called barrel and pincushion distortions produce images that are not high in fidelity compared to the object.

More lens elements = better correction,
but also possibly less light throughput




Figure 4-9

Objective lens designs. Two popular lenses for fluorescence microscopy are shown. Apochromatic lenses may contain 12 or more lens elements to give bright, flat images with excellent color correction across the visual spectrum. Fluorite lenses have fewer lens components and produce sharp, bright images. These lenses exhibit excellent color correction and transmit UV light.


Achromat = corrected for 2 colors

Apochromat = corrected for 3 colors

Plan = flat-field (although not always to full limits of field of view)

TEST BEFORE BUYING!!


Resolution




LABEL-FREE IMAGING

Contrast Enhancing Techniques


Bright
Field



Phase
Contrast




Differential
Interference
Contrast




Phase Contrast & DIC

Phase Contrast Imaging of Transparent Thin Specimens



Differential Interference Contrast Imaging of Transparent Thin Specimens



Differential Interference Contrast Optical Pathways and Components









Figure 3




Quantitative Phase Imaging



Quantitative Phase Tomography





Quantitative Phase Tomography



FLUORESCENCE IMAGING

Fluorescence Basics




<https://www.horiba.com/int/scientific/technologies/fluorescence-spectroscopy/what-is-fluorescence-spectroscopy/>

Fluorescent Probes

Fluorescent proteins




GFP




Organic Dyes

a Shifting spectral properties





b Improving photostability and brightness





Lavis LD. 2017.
Annu. Rev. Biochem. 86:825–43


Anatomy of Fluorescence Microscope



Multi-color Imaging




Multi-color Imaging



6-color LLS imaging of organelles in COS7 cell

ENHANCED RESOLUTION AND SUPER- RESOLUTION IMAGING

Image formation in fluorescence microscope



$$I(\mathbf{r}) = U(\mathbf{r}) \otimes [E(\mathbf{r})c(\mathbf{r})]$$

$I(\mathbf{r})$: *Image*

$U(\mathbf{r}) = PSF(\mathbf{r})$: *Point-spread Function*

$E(\mathbf{r})$: *Excitation Intensity*

$c(\mathbf{r})$: Fluorophore distribution in sample

Wide-Field Imaging

$$E(\mathbf{r}) = \text{const}$$



$$I(\mathbf{r}) = \text{PSF}(\mathbf{r}) \otimes c(\mathbf{r})$$

$$\text{OTF}(\mathbf{k}) = \text{FT}(\text{PSF}(\mathbf{r}))$$

$$\text{OTF}(\mathbf{k}=0) = 0$$

No optical sectioning

OTF support
kmax
NA=nsin(α)



Confocal Microscopy




Image on pinhole plane $\mathbf{s}$, given beam on sample at $\mathbf{r}$


$$I(\mathbf{r}, \mathbf{s}) = \int d\mathbf{r}' U(\mathbf{r} - \mathbf{r}' + \mathbf{s}) [E(\mathbf{r} - \mathbf{r}') c(\mathbf{r}')]]$$

Confocal image, integrate over pinhole (radius = a):


$$I_{conf}(\mathbf{r}) = \int d\mathbf{s} * I(\mathbf{r}, \mathbf{s}) =$$

$$\boxed{\int_{s=0}^a ds U(\mathbf{r} - \mathbf{r}' + \mathbf{s})} \int d\mathbf{r}' [E(\mathbf{r} - \mathbf{r}') c(\mathbf{r}')] =$$

$$PSF_a(\mathbf{r}) \otimes [E(\mathbf{r}) c(\mathbf{r})]$$



Epifluorescence



Confocal

Confocal PSF and OTF


For small pinhole (radius = $a \rightarrow 0$):

$$I_{conf}(\mathbf{r}) = \int d\mathbf{s} * I(\mathbf{r}, \mathbf{s}) =$$

$$\int d\mathbf{r}' U(\mathbf{r} - \mathbf{r}') E(\mathbf{r} - \mathbf{r}') c(\mathbf{r}')$$

$$PSF_{eff}(\mathbf{r}) = PSF_{det}(\mathbf{r}) * PSF_{excit}(\mathbf{r})$$

$$OTF_{eff}(\mathbf{k}) = OTF_{det}(\mathbf{k}) \otimes OTF_{excit}(\mathbf{k})$$




Wide-Field


Confocal

PSF

OTF



Confocal performance vs. pinhole diameter



Axial:

$$FWHM_{tot,axial} = \frac{0.64 \cdot \bar{\lambda}}{(n - \sqrt{n^2 - NA^2})} \quad (7)$$

If $NA < 0.5$, equation (2) can be approximated by:

$$\approx \frac{1.28 \cdot n \cdot \bar{\lambda}}{NA^2} \quad (7a)$$

Lateral:

$$FWHM_{tot,lateral} = 0.37 \frac{\bar{\lambda}}{NA} \quad (8)$$

Axial:

$$FWHM_{ill,axial} = \frac{0.88 \cdot \lambda_{exc}}{(n - \sqrt{n^2 - NA^2})} \quad (2)$$

n = refractive index of immersion liquid
 NA = numerical aperture of the objective
 λ_{exc} = wavelength of the excitation light

If $NA < 0.5$, equation (2) can be approximated by:

$$\approx \frac{1.77 \cdot n \cdot \lambda_{exc}}{NA^2} \quad (2a)$$



Lateral:

$$FWHM_{ill,lateral} = 0.51 \frac{\lambda_{exc}}{NA} \quad (3)$$

Theoretically 1.4x better resolution
 But at expense of signal!

Increase in axial resolution – 4Pi/I5M

2PE 4Pi type C



I5M





Image Scanning and Photon reassignment microscopy

Image on pinhole plane $\mathbf{s}$, given beam on sample at $\mathbf{r}$ (4D data)

$$I(\mathbf{r}, \mathbf{s}) = \int d\mathbf{r}' U(\mathbf{r} - \mathbf{r}' + \mathbf{s}) E(\mathbf{r} - \mathbf{r}') c(\mathbf{r}')$$

$U_{\text{eff}}(\mathbf{r})$, peak is between $\mathbf{r} + \mathbf{s}$ and $\mathbf{r}$.

For photons detected at $\mathbf{s}$, most likely emitter position is at $\mathbf{r} + \mathbf{s}/2$.

Reconstruction:

- shift signal to most probable position
- sum over $\mathbf{s}$

Image Scanning and Photon reassignment microscopy

$$I_{IMS}(\mathbf{r}) = \int d\mathbf{s} * I(\mathbf{r} - \mathbf{s}/2, \mathbf{s}) = \int d\mathbf{r}' \int d\mathbf{s} U(\mathbf{r} - \mathbf{r}' + \mathbf{s}/2) E(\mathbf{r} - \mathbf{r}' - \mathbf{s}/2) c(\mathbf{r}')$$

Substitute: ($s = 4v - 2r + 2r'$)


$$U_{IMS}(\mathbf{r}) = 4 \int d\mathbf{v} U(2\mathbf{v}) E(2\mathbf{r} - 2\mathbf{v})$$

$$\widetilde{U_{IMS}}(\mathbf{k}) = \left[2\tilde{U}\left(\frac{\mathbf{k}}{2}\right) \right] * \left[\tilde{E}\left(\frac{\mathbf{k}}{2}\right) \right]$$

2x resolution

“attenuation factor”

Image Scanning Microscopy






Measure: $I_{IMS}(\mathbf{r})$

Fourier Transform:






$$\widetilde{I_{IMS}}(\mathbf{k})$$

Wiener filter: $\widetilde{I_{IMS}}(\mathbf{k})W(\mathbf{k})$, $W(\mathbf{k}) = \frac{\tilde{U}(\frac{\mathbf{k}}{2})}{\epsilon + \tilde{U}(\frac{\mathbf{k}}{2}) * \tilde{E}(\frac{\mathbf{k}}{2})}$ $\epsilon=0.01$


Reverse FT



Optical Photon Reassignment Microscopy



Structured Illumination Microscopy



Illumination

$$I_{\theta,\varphi}(r) = \frac{I_0}{2} [1 + \cos(2\pi k_\theta \cdot r + \varphi)]$$


Image

$$D(r) = (s(r) \cdot I_{\theta,\varphi}(r)) \otimes h(r)$$

$$D_{\theta,\varphi}(k) = [S(k) \otimes I_{\theta,\varphi}(k)] \cdot H(k)$$

$$= \frac{I_0}{2} [S(k) \cdot H(k) + \underbrace{\frac{1}{2} e^{-i\varphi} S(k - k_\theta) \cdot H(k)}_{\text{Information previously inaccessible}} + \underbrace{\frac{1}{2} e^{i\varphi} S(k + k_\theta) \cdot H(k)}_{\text{Information previously inaccessible}}]$$

Information previously inaccessible



Images with
additional
information

Individual Information
components

Final reconstruction

Acquire 3 images ($\phi=0, 2\pi/3, 4\pi/3$), solve linear system

$$D_{\theta,\varphi_1}(k) = \frac{I_0}{2} [S(k) \cdot H(k) + \frac{1}{2} e^{-i\varphi_1} S(k - k_\theta) \cdot H(k) + \frac{1}{2} e^{i\varphi_1} S(k + k_\theta) \cdot H(k)]$$

$$D_{\theta,\varphi_2}(k) = \frac{I_0}{2} [S(k) \cdot H(k) + \frac{1}{2} e^{-i\varphi_2} S(k - k_\theta) \cdot H(k) + \frac{1}{2} e^{i\varphi_2} S(k + k_\theta) \cdot H(k)]$$

$$D_{\theta,\varphi_3}(k) = \frac{I_0}{2} [S(k) \cdot H(k) + \frac{1}{2} e^{-i\varphi_3} S(k - k_\theta) \cdot H(k) + \frac{1}{2} e^{i\varphi_3} S(k + k_\theta) \cdot H(k)]$$



$$\begin{bmatrix} S(k) \cdot H(k) \\ S(k - k_\theta) \cdot H(k) \\ S(k + k_\theta) \cdot H(k) \end{bmatrix} = \frac{2}{I_0} M^{-1} \begin{bmatrix} D_{\theta,\varphi_1}(k) \\ D_{\theta,\varphi_2}(k) \\ D_{\theta,\varphi_3}(k) \end{bmatrix} \quad M = \begin{bmatrix} \frac{1}{2} e^{-i\varphi_1} & \frac{1}{2} e^{i\varphi_1} \\ \frac{1}{2} e^{-i\varphi_2} & \frac{1}{2} e^{i\varphi_2} \\ \frac{1}{2} e^{-i\varphi_3} & \frac{1}{2} e^{i\varphi_3} \end{bmatrix}$$

Reconstruct by shifting components in Fourier Space


$$\widehat{S}(k) = \frac{\sum_{\theta,m} H_m^*(k + mk_\theta) \widehat{D}_{\theta,m}(k + mk_\theta)}{\sum_{\theta',m'} |H_{m'}(k + m'k_{\theta'})|^2 + w^2} A(k)$$

$$\widehat{D}_{\theta,m}(k + mk_\theta) = S(k + mk_\theta) \cdot H(k), \quad m=-1,0,1$$

3D Structured Illumination Microscopy (SIM)






Structured Illumination Microscopy




Caveolae dynamics in a COS-7 cell

I5S (SIM + 4Pi)




Stimulated emission depletion (STED) microscopy








Resolution scaling in STED microscopy



$$d \approx \frac{\lambda}{2NA\sqrt{1+\zeta}}, \zeta = I_{STED}/I_{STED}^{Sat}$$



Combination of STED with 4Pi interferometry




Single-molecule localization microscopy (SMLM)




Discoveries by SMLM


a Nucleosome clutches



b Nuclear pore complexes



c Focal adhesion nanoarchitecture



d Actin rings





Image-based Single-molecule Localization




Fundamental limit: **signal-to-noise** (“uncertainty principle”)


Theory: Thompson et al (2002); Molecular motor tracking: Yildiz et al (2003).

3D Localization approaches


Bi-Plane






PSF Engineering




Interferometry




Advanced Techniques: Modulation-enhanced Localization



~ 1 nm 3D localization precision (1σ)






Modulation enhanced localization





IMAGING LIVE CELL AND MULTI-CELLULAR SYSTEMS


Multi-photon Microscopy





Multi-photon Microscopy





Trade-offs in imaging live biology



Selective-plane Illumination Microscopy (SPIM)




Gaussian Beam propagation





$$z_R = \frac{\pi w_0^2 n}{\lambda}$$

Trade-off between spatial resolution and FOV





SPIM with “Non-diffracting” Bessel beams




Problem: significant side-lobes and out-of-focus excitation


Bessel Beam Arrays and Bounded Optical Lattices





Bounded 2D Optical Lattices







SPIM with Bounded Optical Lattices (Lattice Light-Sheet Imaging)



Lattice Light-Sheet Imaging




3D Interferometric LLS Imaging



Increased 3D resolution and OTF support

Selective-plane Illumination Microscopy (SPIM)



Selective-plane Illumination Microscopy (SPIM)







Image Degradation in thick samples




Source: Flickr




Closed-loop AO with direct wavefront sensing





Shack-Hartmann wavefront sensor




Adaptive Optics Imaging



Adaptive Optics Imaging (Wide-field)



Adaptive Optics Imaging



Adaptive Light-sheet Imaging

